

# Conventional Macintosh Blade

Crimson WARMLITE & LEDLITE series are designed with quick to replace lamp at the front of blades. These blades has low profile with cost effective features and are simply autoclavable for more than 2000 times at specific temperature. This contributes mostly to the elimination of cross-infection. Crimson WARMLITE & LEDLITE comes with German quality bulbs which delivers matchless illumination and bright shine intensity. The bulbs are available in two lights LED & Halogen. All blades are in accordance of ISO 7376/1.



- With reliable 2.5V light source for good illumination of the examination area.
- Autoclavable upto 134°C approximately 2000 times.
- No adverse effect in using ETO, steam autoclaving.
- All blades are made from Austenitic Stainless steel thereby making them rust proof.
- Specially designed blades for easy intubation.
- Round edges prevent from injury risk.
- Compatible with ISO 7376/1 standard handles.
- Matt finished to reduced glare.
- Having good light transmission of 3500 lux with 2.5V handle.
- Can sustain plasma sterilization. (STERRAD 100NX)

**SIZE**

CI-1201-00	00
CI-1201-0	0
CI-1201-1	1
CI-1201-2	2
CI-1201-3	3
CI-1201-4	4
CI-1201-5	5